

XLIFE

Step 1: Download the APP

iOS: Search "Xlife" in App Store Android: Search "Xlife" in Google Play

Or scan the QR code to download the app



Step 2: Press and hold the power button to turn on the ear cleaning rod



Step 3: Connect product WiFi



Connect Product Wi-Fi "Xlife-XO-XXXXXX"

Back to Xlife app Enter and start to use it



Ear Otoscope Charging Cable User Manual x 1 x 1

Accessory



Silicone Ear Spoon Cover x 8

Note: Do not use the charging equipment with output voltage exceeding 5V

Troubleshooting

| Problem | Solution | | | |
|---|---|--|--|--|
| No response after startup | Please use the incidental charging cable to charge. The charging time is 30mins | | | |
| The app exits quickly | Delete the existing app and scan QR to reinstall it | | | |
| Indicator light flashes | Please use the incidental charging cable to charge. The charging time is 30mins | | | |
| Abnormal display in app | Interrupt the connection, restart the product to connect again | | | |
| If you still have problems with the above procedures, | | | | |

please contact our Customer Service Representatives

directly for further assistance

Hazardous substance

| Part name | Hazardous material | | | | | | |
|--------------------------------|--------------------|------|------|----------|-------|--------|--|
| | (Pb) | (Hg) | (Cd) | (Cr(VI)) | (PBB) | (PBDE) | |
| Printed Circuit Components | × | 0 | 0 | 0 | 0 | 0 | |
| Case | 0 | 0 | 0 | 0 | 0 | 0 | |
| Plastic rubber accessory | 0 | 0 | 0 | 0 | 0 | 0 | |

This form is formulated in accordance with SJ/T 11364

- ○: represent the hazardous materrial ratio is below GB/T26572-2011 x: represent the hazardous materrial ratio is over GB/T26572-2011

Product parameters

Product Name

Pixels : 3 Mega-pixels

: 100mAh

: 2.4Ghz

Depth of field : 10-50mm : 3.5mm Lens diameter

Image transmission rate : 20fps

Battery life : 30 minutes

Charging time : 30 minutes

: 5V 0.5A Input current : 3**-**axis

Network standard : IEEE 802.11b/g

Working temperature : -10- 45 C

: 11g Rod weight

: CMOS Image sensor

Attention

Working frequency

Battery

Gyroscope

- 1. When cleaning lens, please wipe it carefully with professional alcohol cotton swabs.
- 2. Please observe your surroundings before use. Stay at least 2ft away from people, animals, or moving objects when using. Don't use in a moving vehicle or other unstable environments. Please don't move around when using.
- 3. This product is unsuitable for children under 3 years old.
- 4. This product is not suitable for children to use alone.
- 5. Do not place in a humid environment for a long time or
- 6. There is a built-in rechargeable lithium-ion battery. Even if you do not use it, you should charge it once a month to avoid long-term battery damage.
- 7. Please use a 5V charger to charge the device. Do not use this product while being charged.
- 8. Avoid direct exposure to the sun. 9, Do not put heavy objects on the product.
- 10. It is normal for the temperature to rise slightly when using (up to 95°F).



www.portronics.com For Support

Visit: www.portronics.com **Email:** help@portronics.com

Contact No.: +91 9555245245, 011-42413131 (Working Hrs.: Monday - Saturday 10:00am - 6:00pm)



For E-waste & ROHS compliance refer to http://www.portronics.com/ewaste.html

EAR CLEANING OTOSCOPE